★★★
TRISTAR

EN | Instruction manual

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PT | Manual de utilizador

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SV | Bruksanvisning

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SK | Návod na použitie

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EN Instruction manual

SAFETY

- By ignoring the safety instructions the manufacturer can not be hold responsible for the damage.
- The appliance must be placed on a stable, level surface.
- The body fat analysis scale is not suitable for the following persons:

- Children under 10 years and adults over 85 years
- Persons with symptoms of fever, oedema, or osteoporosis
- Persons undergoing treatment for dialysis
- Persons who use a pacemaker
- Persons who take cardiovascular medication
- Pregnant women
- Sportspersons who perform over 5 hours of intensive sporting activities each day
- Persons with a heartbeat lower than 60 beats/minute

Batteries

- Do not expose the battery to high temperatures or direct sunlight.
- Never throw batteries into the fire. There is a danger of explosion!
- Keep batteries away from children. They are not a toy!
- Do not open the batteries by force.
- Avoid contact with metallic objects.(Rings, nails, screws et cetera)there is a danger of short-circuiting!
- As a result of a short-circuit batteries may heat up consider-ably or even catch fire. This may result in burns.
- For your safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks, do not rub your eyes with your hands. If in contact with the liquid, the hands should be washed thoroughly and eyes rinsed out with clear water, if symptoms persist please consult a doctor.

PARTS DESCRIPTION

- Contact surface
- Display
- Down button
- Up button
- Set button

Display symbol	Meaning
Kg./lb/st	Body weight
P0-9	User memory space
	Male
	Female
CM(FT;in)	Body height
Age	Age
FAT	Body fat display
TBW	Body water display
MUS	Muscle mass display
BONE	Bone mass display
%	Percentage body fat / body water / muscle mass/ bone weight

BEFORE THE FIRST USE

- Take the appliance and accessories out the box. Remove the stickers, protective foil or plastic from the device.
- The device can only function correctly if you stand on the contact surface with bare, dry feet!
- The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body.
- The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor we measure is the relationship between body weight and associated body water percentages as well as on other biological data(age, gender, body weight).
- There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values, fat is a factor that puts serious strain on the body, in particular the cardio-vascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.
- This signal is generally safe! However, please read our special safety instructions!

USE

- The scales are factory set with the units"kg" and "cm". To change the settings to "lb" or "kg" and "Ft.in", proceed as follows: On the back side of the scale beside the battery compartment a switch can be found, push the below switch, choose between the measuring units kg/lb/st.

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

Saving Personal Date

- The scale has default values that you can change, as described below.
- Carry out the following steps in sequence and within a few seconds, otherwise the device will automatically cancel the input mode.
- Press the button SET. The memory space will flash on the display (P0). Confirm the setting by pressing SET or select another memory space by pressing the UP/DOWN button and then SET to go to the next setting.
- The symbol for"male" or "female" will flash. Confirm the setting by pressing SET or use the UP/DOWN button to choose the symbol for your gender and then press SET to go the next setting.
- The height setting will flash. Select your height by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The setting AGE will flash. Select your age by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.

Weigh on based on personal settings

- The device can only function correctly if you stand on the device with bare,dry feet. If necessary wipe the contact surface to improve the results.
- Due to the significant fluctuation of body fat the optimal testing time is between 6 p.m. And 8 p.m, since the body fat remains very consistant during this period of time
- Press the SET button, the user number will flash (P0-P9). Use the UP/ DOWN button to select your preset data.
- Wait for the user number to stop flashing.
- Stand on the contact surfaces with both feet.
- Firstly your weight will flash on the display.
- After a few seconds, the following results will be shown:
 - the body fat content in% (FAT symbol)
 - the results analysis in the form of a bar
 - the body water content in % (TBW symbol)
 - the muscle mass content in % (MUS symbol)
 - the bone mass content in %(BON symbol)
- Error
- If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.
- 8888 = Self test (when switching on) to check if the numeric display functions properly.
- = The scale is not placed properly.
- Err = The Percentage of body fat content lies above or below the limits of the measuring scale
- EEEE = The weight lies above the allowed maximum weight.

CLEANING AND MAINTENANCE

- This product is guaranteed for 24 months granted. Your warranty is valid if the product is used in accordance to the instructions and for the purpose for which it was created. In addition, the original purchase (invoice, sales slip or receipt) is to be submitted with the date of purchase, the name of the retailer and the item number of the product.
- For the detailed warranty conditions, please refer to our service website: www.service.tristar.eu

GUARANTEE

- This product is guaranteed for 24 months granted. Your warranty is valid if the product is used in accordance to the instructions and for the purpose for which it was created. In addition, the original purchase (invoice, sales slip or receipt) is to be submitted with the date of purchase, the name of the retailer and the item number of the product.
- For the detailed warranty conditions, please refer to our service website: www.service.tristar.eu

ENVIRONMENT

- This appliance should not be put into the domestic garbage at the end of its durability, but must be offered at a central point for the recycling of electric and electronic domestic appliances. This symbol on the appliance, instruction manual and packaging puts your attention to this important issue. The materials used in this appliance can be recycled. By recycling of used domestic appliances you contribute an important push to the protection of our environment. Ask your local authorities for information regarding the point of recollection.

Support

You can find all available information and spare parts at service.tristar.eu!

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